

beyontics

Peakprovider

Defining Unknown Peaks in
Empower®

beyontics **PeakProvider 2.0 for Empower®**

A smart way of analyzing unknown peaks in Empower®

beyontics PeakProvider assists in analyzing unknown peaks by accelerating the evaluation of chromatograms with dozens of unidentified components and making results more reliable. Transfer to external calculation programs will belong to the past.

Analyzing unknown peaks in Empower. Do you recognize yourself?

When having lots of unknown peaks in your data, difficulties might occur for defining suitable peak tables in your processing method. Without these entries inter-sample calculations are not possible. Creating the best matching peak entries takes time and might cause difficulties. A manual check for each chromatogram within your data set against once defined peak entries needs to be carried out repeatedly, retention times need to be adjusted. Getting the best matching definitions in your table requires substantial time and effort. Just having a few unknowns might not be a big issue, but what about having dozens of them?

An automated solution. beyontics PeakProvider.

The PeakProvider is a smart and easy-to-use tool with a huge effect when lots of unknown peaks are present in your data on which inter-sample calculation needs to be carried out. By the application of a subtle algorithm, the tool determines identification criteria that matches all the chromatgrams within your data set, not just a single chromatogram. A unique name is assigned, by which each component can be identified clearly from injection to injection. No manual adjustment of your retention times is required anymore.

Easily integrated into Empower. How does it work?

With the PeakProvider beyontics provides a tool that is integrated in the Empower software and is

especially useful for analyzing sample solutions containing a high number of unknown peaks. Additional peak names are added to the component table of an already existing processing method in an automated way, based on a user-defined selection of chromatograms.

Easy to use. The PeakProvider workflow.


Peak identification is performed automatically by using either the absolute or relative retention time as part of the assigned peak name. Any prefixes or suffixes can be configured by the user.

Step 1:

Process your sample set in the usual way for generating a result set.

Step 2:

Apply the PeakProvider on the result set and define the kind of the desired identification for unknown peaks.



Step3:


A new result set is reprocessed by the PeakProvider. The applied processing method now contains peak entries for all identified unknown peaks in the component table.

Integration		Smoothing/Offset		Purity	PDA Library Search	Components	Impurity	Peak
Average By	None						Update RT	Never
RT Window (%)	5.00	CCalRef1		PP208_01				
<input checked="" type="checkbox"/> Include Internal Std Amounts in % Amount Calculation								
Sample Value Type	Amount	Auto Peak Label		RT Reference Used to Unnamed Peaks by Rt				

Name	Component Type	Peak Label	Retention Time (min)	RT Window (min)	Peak Match	Y Value	X Value
1 OQ # 0,224			1.532	0.077	Closest	Area	Amount
2 OQ # 0,251			1.725	0.086	Closest	Area	Amount
3 OQ # 0,364			2.499	0.125	Closest	Area	Amount
4 OQ # 0,377			2.609	0.130	Closest	Area	Amount
5 OQ # 0,925			6.354	0.318	Closest	Area	Amount
6 PP208_01			6.900	0.345	Closest	Area	Amount
7 OQ # 1,432			9.901	0.495	Closest	Area	Amount


Names for all unknown peaks are added to the processing method. The name contains for example the relative retention time, that was calculated as a mean value based on all chromatograms that have been considered during calculation. Thus, the PeakProvider makes sure, that the same component can be discovered in every chromatogram of the result set.

If the user has decided, that a new result set should be created by applying the updated processing method, all peaks in the reprocessed result set are clearly identified by a peak name. Summary evaluation in report tables or by custom fields are now possible for all known and unknown peaks.



The functions of the Peakprovider at a glance:

- » *integrated into Empower as a toolkit application*
- » *intelligent algorithm to identify unknown peaks and assign unique peak names to the same components by choosing either relative or absolute retention time values*
- » *automated detection of the most suitable section width to clearly identify same components of the operated chromatograms and assign unique component names*
- » *a naming convention for the unknown peaks can be configured by defining prefixes or suffixes*
- » *flagging of components added by the tool is possible*
- » *decision whether or not flagged component names should be removed from the component table prior to a new identification step*
- » *beside the retention time any customized peak field can be defined as a value, on which the identification should be based.*
- » *full support of the Empower audit trail function*



Still not convinced? Get your free trial version now!

We offer a 30 days trial version for free. Please don't hesitate to contact us. For contact details, please look below.

beyonitics Headquarter

beyonitics GmbH

Altonaer Straße 79-81

D- 13581 Berlin

Phone: +49 30 36432319

Fax: +49 30 36431977

E-Mail: info@beyonitics.com

Registry Court:

Berlin- Charlottenburg, HRB 74211

VAT-ID: DE185286048

beyonitics Regional Organisations and Agencies

Germany and EU

beyonitics Europe GmbH

Phone: +49 30 36432319

E-Mail: info@beyonitics.com

North and South America

beyonitics America Inc.

Phone: +1 857 2142098

E-Mail: america@beyonitics.com

Asia-Pacific

beyonitics Asia-Pacific

Phone: +63 939 8986264

E-Mail: asia-pacific@beyonitics.com

Switzerland and other regions

beyonitics Schweiz GmbH

Phone: +41 61 2254355

E-Mail: schweiz@beyonitics.com